



ADD050

45 ...
.052 (...)



... #50 ...

Piece #3

Dr. Boren
5/11/71

with needle + shield } not only
with shield }
screw }
(stated)



rubber 1/32 thickness

stainless steel .025 thick




8/32 TME


make 6 each

8-32 thread

* Same has been .750 for this dimension.



opening #18 gauge
ID .040



MAKE-2

Perfusion piece #3

(1 Prototype)

- 6 - Rubber seals
- 6 - 9/16 pressure disk
- 6 - 8/32 slotted screws

Perfusion piece #4